

K30 Gastight Syringe




- Standard 1/4-28UNF thread, good tightness
- Wetted Material: Borosilicate(glass), PTFE(piston /syringe outlet) , high precision, corrosion resistance
- Maximum pressure rating 0.68Mpa(Bar)
- Easy clean and sterilization, reusable
- Service Life: 3 million times no leakage (water)
- Flexible connections to meet different applications

RZ-K30 - ^①5ML - W - 1 - U
 Model No. Volume Male thread Single channel Thread 1/4-28UNF

5ML	250μl
2.5ML	125μl
1.25ML	50μl
500μl	25μl

Dimension (Unit: mm)




5ml




500μl




250μl



125μl




50μl



2.5ml

Other specifications share the same dimension with 5 ml



1.25ml

Other specifications share the same dimension with 5 ml



25μl

Other specifications share the same dimension with 5 ml